



GE Healthcare

Panda™ Warmer

For when baby can't be with mom

An integrated approach to confident care

The best place for your newborns is in the loving arms of their moms and dads. Sometimes, though, they need extra care and attention—the kind of care and attention only you can provide.

You want to give every baby, no matter how small the size or big the challenge, the chance for the best start at life. The most important tools you have is your clinical know-how paired with your empathy to give newborns the chance they deserve.

Panda Warmers from GE Healthcare give you the all-in-one solution you need to help babies go home healthy. Panda Warmers are developed with an emphasis on patient-centered care to create a more inclusive and collaborative environment for care providers and their patients. In addition, Panda Warmers provide the integrated technology solutions to help you deliver responsive, appropriate care.



One integrated solution

Labor and delivery providers are masters of juggling multiple tasks, while staying vigilant for any number of possible challenges. The features of the Panda Warmers are designed to allow providers to focus on their most important task—caring for newborn patients and their moms.



Simple

A bright idea

The optional integrated procedure light is simple to position and lets you easily focus the light where it is needed. The standard observation light allows you to adjust the light intensity, so that you can do the work you need to do without creating any unnecessary discomfort for the newborn.

Brilliant, full-color controls

The control panel & full-color display are accessible and easy to use. Patient information is vividly displayed on the large screen and is easy to view.

Dovetail dual rails

The GE Healthcare vertical dovetail dual rail system is incorporated for easy addition of accessories and ancillary equipment to simplify your equipment management.



Efficient

Warm-up with no nuisance alarms

When turning on the Panda Warmer, the automatic warm-up mode maximizes warming to help ensure the bed is warm when the newborn arrives. By eliminating nuisance alarms, we help keep the environment quiet and stress-free for caregivers and parents.

Heat where you need it

GE Healthcare's innovative Hourglass Heat Profile bathes the baby in soothing warmth, with uniform coverage across the mattress, while helping to keep the caregiver cool and comfortable.



Intelligent

Hands-free alarm silence

Maintain a clean environment and keep the area quiet and soothing by simply moving your hand in front of the sensor to silence the alarm.

Accurate temperature control

Powered by the same closed-bed, heat-controlled algorithm as our trusted Giraffe™ OmniBed™ Carestation.

Next level of integration

Make the most of your investment and provide the best possible care with these optional features for your Panda Warmer.

Integrated resuscitation

Available with minimal setup to meet NRP and ILCOR guidelines, and helps to standardize resuscitation protocols across the perinatal care area. Integrated gas cylinder holders and tank pressure gauges easily visible to enable smooth workflow and transport convenience. Integrated t-piece system with a built in manometer and integrated air/oxygen blender so all your tools are at your fingertips.

Integrated SpO₂ and ECG Heart rate

The integrated ResusView™ Solution, available with Panda iRes warmers, provides accurate ECG heart rate and SpO₂ blood oxygen levels within seconds, allowing you to provide timely and appropriate interventions and the best care possible for babies and families.

Integrated scale

Baby's first weight is available right at the bedside. There is no need to move the baby to another surface and elevates comfort with a warm surface rather than a cold scale. So cute, parents might even take a picture!

X-ray tray

With an integrated x-ray cassette tray, there's no need to move the baby; the Panda Warmer maintains uninterrupted thermal support during mobile radiography exams.

Uninterruptible Power Supply (UPS)*

Quality of care is never compromised, as the UPS battery keeps the baby safe and secure from power interruptions for 15 minutes.**



Three options. Infinite potential.

No matter which Panda Warmer is best for your labor and delivery facility, you can be confident that each one is designed to empower you to give your newborns the care they need, when they need it most.



Panda Wall Mount Warmer

The Panda Wall Mount Warmer provides a wide range of features while complementing the design and architecture of your labor and delivery room. Perfect for enclosure in custom cabinetry, the Panda Wall Mount Warmer can be out of sight when not needed and instantly available when it is needed.



Panda Freestanding Warmer

The Panda Freestanding Warmer works with your own bed or bassinet and can be easily moved throughout your facility, improving your staff workflow and minimizing disruptions to your tiniest patients and families. The Panda Freestanding Warmer offers a complete range of features combined with exceptional versatility.



Panda iRes Warmer

The Panda iRes Warmer provides an integrated, transportable solution to help you provide the strongest start for newborns in your care. With ResusView technology, the Panda iRes Warmer can be configured with full resuscitation capabilities, ECG heart rate and SpO₂ monitoring, meeting clinical standards in a single configuration.

At the heart of the golden hour

You give the world's newborns a strong start at life, and we're here to support you with reliable, advanced technology that's also incredibly easy to use. GE Healthcare's Panda Warmer ResusView solution, available exclusively with the Panda iRes Warmer, helps give your staff clinical confidence in those vital moments.



Heart rate in 6 seconds***

An accurate heart rate assessment immediately after birth is critical to the correct guidance of resuscitation efforts. With ResusView, ECG heart rate and blood oxygen saturation level (SpO₂) are integrated into your warmer.

ResusView also provides a convenient view of this data with a continuous ECG heart rate number and waveform displayed on the Panda Warmer screen, something palpation or auscultation won't provide. Integrated, the tools are in one spot, so they're there when you need them. ResusView heart rate is available as an upgrade to your existing Panda Warmer. Quickly get an accurate heart rate and oxygenation data to help make informed decisions quickly.



Integrated



Turn-key



**Clinical
confidence**





Disclaimers

* Not available in all markets

** Dependent upon load conditions

*** 60-119 bpm

Imagination at work

© 2020 General Electric Company – All rights reserved. GE Healthcare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Healthcare representative for the most current information. GE Medical Systems, Inc., doing business as GE Healthcare. GE Healthcare, a division of General Electric Company. GE, GE Monogram, Panda, and ResusView are trademarks of General Electric Company.

www.gehealthcare.com

JB61441XX